



# Pipette / Dispenser



NEXTY series



Peggy



BlackBird



# NEXTY

## Are you stressed out about pipette volume setting??

NEXTY has 2 kind dials of volume setting. These dials relieve your stress from volume setting because setting volume will be faster.


Mushroom-shaped Push Button

The mushroom-shaped push button reduces the workload of dispensing. You only need very low operating forces because this button has new spring mechanism. This push button functions as dial of minimum volume setting too.

The world's first!!

Triple-Speed Turbo Dial

NEXTY's triple-speed turbo dial performs easy-to-use volume setting. One rotation of Triple Dial is equal to 3.5 rotations of the ordinary push-button dial which is useful when making a large sample volume change.



### Lock lever

Push button and turbo dial's rotation can be locked by this lock lever. Can prevent the deviations in volume setting.



### Eject button

This eject button reduces the workload of tip ejecting.



### Color-coded

NEXTY's push button has different color by each maximum volume.



### Digital indicator

Scale is easy to look at because display is three digit. Can set 1/1000 volume of maximum.

Single Channel

Multi Channel

Autoclavable





*NEXTY Autoclavable Pipettes*

Fully autoclavable (121°C for 20 minutes) !!

You can perform autoclave sterilization, therefore you can use it kept clean without worry about contamination.



- Lock lever
- Color-coded
- Digital indicator



**Tool for Assembly and Disassembly**

Attach the convenient tool for assembly and disassembly when you need maintenance of NEXTY.  
\* Only NEXTY-A5000 is excluded.



**Caution:**  
Do Not Touch immediately after autoclave, because of Hot surface of autoclavable NEXTY.

Cat. No.	Item	Minimal scale	Volume	Accuracy	Precision	Compatible tip	Packaging
NT-A2	NEXTY-A2 Autoclavable Pipette 0.2~2 μL	0.001 μL	0.2 μL	≤ ± 15.0%	≤ 8.0%	<b>204</b> / 207	1 unit
			2 μL	≤ ± 2.5%	≤ 1.0%		
NT-A10	NEXTY-A10 Autoclavable Pipette 1~10 μL	0.01 μL	1 μL	≤ ± 5.0%	≤ 2.5%	204 / <b>207</b>	1 unit
			10 μL	≤ ± 1.0%	≤ 0.5%		
NT-A20	NEXTY-A20 Autoclavable Pipette 2~20 μL	0.01 μL	2 μL	≤ ± 6.0%	≤ 3.0%	<b>705</b> / 703 / 503	1 unit
			20 μL	≤ ± 1.0%	≤ 0.5%		
NT-A100	NEXTY-A100 Autoclavable Pipette 10~100 μL	0.1 μL	10 μL	≤ ± 2.0%	≤ 1.0%	<b>705</b> / 703 / 503 / 505 / 801	1 unit
			100 μL	≤ ± 0.8%	≤ 0.3%		
NT-A200	NEXTY-A200 Autoclavable Pipette 20~200 μL	0.1 μL	20 μL	≤ ± 2.0%	≤ 1.0%	<b>705</b> / 703 / 503 / 505 / 801	1 unit
			200 μL	≤ ± 0.8%	≤ 0.3%		
NT-A1000	NEXTY-A1000 Autoclavable Pipette 100~1000 μL	1 μL	100 μL	≤ ± 2.0%	≤ 0.6%	502 / <b>706</b> / 804 / 806 / 805	1 unit
			1000 μL	≤ ± 0.8%	≤ 0.3%		
NT-A5000	NEXTY-A5000 Autoclavable Pipette 500~5000 μL	5 μL	500 μL	≤ ± 2.0%	≤ 0.6%	<b>401</b>	1 unit
			5000 μL	≤ ± 0.8%	≤ 0.3%		

Note: Each pipette is calibrated with a tip in bold letter shown above. Other tips if used may cause a slight deviation in performance.

### NEXTY Single Channel Pipettes

NEXTY standard models have high operating efficiency.



- Lock lever
- Color-coded
- Digital indicator




Bottom Parts



#### Autoclavable bottom parts

The nose-cone and eject cone can be taken off and autoclaved.

\*O-rings and seal rings can be left in nose-cone - they are autoclavable. The Nose cone can be autoclaved as well. Nose cone can be autoclaved with together them.



Cat. No.	Item	Minimal scale	Volume	Accuracy	Precision	Compatible tip	Packaging
NT-2	NEXTY-2 Single Channel pipette 0.2 ~ 2 $\mu$ L	0.001 $\mu$ L	0.2 $\mu$ L	$\leq \pm 15.0\%$	$\leq 8.0\%$	<b>204</b> / 207	1 unit
			2 $\mu$ L	$\leq \pm 2.5\%$	$\leq 1.0\%$		
NT-10	NEXTY-10 Single Channel pipette 1 ~ 10 $\mu$ L	0.01 $\mu$ L	1 $\mu$ L	$\leq \pm 5.0\%$	$\leq 2.5\%$	204 / <b>207</b>	1 unit
			10 $\mu$ L	$\leq \pm 1.0\%$	$\leq 0.5\%$		
NT-20	NEXTY-20 Single Channel pipette 2 ~ 20 $\mu$ L	0.01 $\mu$ L	2 $\mu$ L	$\leq \pm 6.0\%$	$\leq 3.0\%$	<b>705</b> / 703 / 503	1 unit
			20 $\mu$ L	$\leq \pm 1.0\%$	$\leq 0.5\%$		
NT-100	NEXTY-100 Single Channel pipette 10 ~ 100 $\mu$ L	0.1 $\mu$ L	10 $\mu$ L	$\leq \pm 2.0\%$	$\leq 1.0\%$	<b>705</b> / 703 / 503 / 505 / 801	1 unit
			100 $\mu$ L	$\leq \pm 0.8\%$	$\leq 0.3\%$		
NT-200	NEXTY-200 Single Channel pipette 20 ~ 200 $\mu$ L	0.1 $\mu$ L	20 $\mu$ L	$\leq \pm 2.0\%$	$\leq 1.0\%$	<b>705</b> / 703 / 503 / 505 / 801	1 unit
			200 $\mu$ L	$\leq \pm 0.8\%$	$\leq 0.3\%$		
NT-1000	NEXTY-1000 Single Channel pipette 100 ~ 1000 $\mu$ L	1 $\mu$ L	100 $\mu$ L	$\leq \pm 2.0\%$	$\leq 0.6\%$	502 / <b>706</b> / 804 / 806 / 805	1 unit
			1000 $\mu$ L	$\leq \pm 0.8\%$	$\leq 0.3\%$		
NT-5000	NEXTY-5000 Single Channel pipette 500 ~ 5000 $\mu$ L	5 $\mu$ L	500 $\mu$ L	$\leq \pm 2.0\%$	$\leq 0.6\%$	<b>401</b>	1 unit
			5000 $\mu$ L	$\leq \pm 0.8\%$	$\leq 0.3\%$		


Note: Each pipette is calibrated with a tip in bold letter shown above. Other tips if used may cause a slight deviation in performance.

### NEXTY Single Channel Pipettes LONG NOZZLE

You can use easily for long test tubes!!  
Useful long nozzle type.



- Lock lever
- Color-coated
- Digital indicator
- Autoclavable bottom parts



NEXTY Long Nozzle / Multi Channel Pipettes

Cat. No.	Item	Minimal scale	Volume	Accuracy	Precision	Compatible tip	Packaging
NTLG-200	NEXTY Single Channel pipette LONG NOZZLE 20 ~ 200 $\mu$ L	0.1 $\mu$ L	20 $\mu$ L	$\leq \pm 2.0\%$	$\leq 1.0\%$	<b>705</b> / 703 / 503 / 505 / 801	1 unit
			200 $\mu$ L	$\leq \pm 0.8\%$	$\leq 0.3\%$		
NTLG-1000	NEXTY Single Channel pipette LONG NOZZLE 100 ~ 1000 $\mu$ L	1 $\mu$ L	100 $\mu$ L	$\leq \pm 2.0\%$	$\leq 0.6\%$	502 / <b>706</b> / 804 / 806 / 805	1 unit
			1000 $\mu$ L	$\leq \pm 0.8\%$	$\leq 0.3\%$		


Note: Each pipette is calibrated with a tip in bold letter shown above. Other tips if used may cause a slight deviation in performance.

### NEXTY Multi Channel Pipettes



- Lock lever
- Color-coated
- Digital indicator

#### Eight-channel multi-suspension system



Eight lanes are respectively equipped with independent suspension to ensure accurate tip insertion.



Cat. No.	Item	Minimal scale	Volume	Accuracy	Precision	Compatible tip	Packaging
NT-810	NEXTY 8 Multi Channel pipette 1 ~ 10 $\mu$ L	0.01 $\mu$ L	1 $\mu$ L	$\leq \pm 8.0\%$	$\leq 5.0\%$	<b>207</b> / 204	1 unit
			10 $\mu$ L	$\leq \pm 2.0\%$	$\leq 1.0\%$		
NT-8200	NEXTY 8 Multi Channel pipette 20 ~ 200 $\mu$ L	0.1 $\mu$ L	20 $\mu$ L	$\leq \pm 8.0\%$	$\leq 1.0\%$	<b>705</b> / 703 / 503 / 505 / 801	1 unit
			200 $\mu$ L	$\leq \pm 1.0\%$	$\leq 0.5\%$		

Note: Each pipette is calibrated with a tip in bold letter shown above. Other tips if used may cause a slight deviation in performance.

### NEXTY Fixed Volume Pipettes

The sample volume can be customized upon request.

NEXTY Fixed Volume Pipettes



- Color-coated
- Autoclavable bottom parts



#### Eject button

The unique design of our eject button provides easy and comfortable ejecting and dispensing.




#### Smooth & Steady stroke

Our NEXTY series of pipette passing half a million repetitive stroke tests, always have a smooth and precise stroke by adopting stable spring system and PBT internal plunger which is highly resistant to chemicals.

#### Volume Change Tool Kit

for NEXTY Fixed Volume Pipette

The pipette can be customized as an “original fixed pipette” by setting the desired capacity.



A video showing the procedure is available on our website.

[watsonbiolab vol-change](#)

Cat. No.	Item	Unit
NT-FVC	Volume Change Tool Kit	1 unit

Cat. No.	Item	Volume	Accuracy	Precision	Compatible tip	Packaging
NT-F10	NEXTY-F10 fixed volume pipette 10 µL	10 µL	≤ ± 1.0%	≤ 0.5%	204 / 207	1 unit
NT-F20	NEXTY-F20 fixed volume pipette 20 µL	20 µL	≤ ± 1.0%	≤ 0.5%	703 / 705	1 unit
NT-F100	NEXTY-F100 fixed volume pipette 100 µL	100 µL	≤ ± 0.8%	≤ 0.3%	703 / 705	1 unit
NT-F200	NEXTY-F200 fixed volume pipette 200 µL	200 µL	≤ ± 0.8%	≤ 0.3%	703 / 705	1 unit
NT-F1000	NEXTY-F1000 fixed volume pipette 1000 µL	1000 µL	≤ ± 0.8%	≤ 0.3%	706 / 804	1 unit
NT-F5000	NEXTY-F5000 fixed volume pipette 5000 µL	5000 µL	≤ ± 0.8%	≤ 0.3%	401	1 unit





NEXTY Pipette stand

**NEXTY Pipette stand**

**Pipette Carrousel Stand for NEXTY series.**

Pipette stand combining functionality and excellent design.

This pipette stand is compatible with some pipettes from other manufacturers.



\* The product package does not include pipettes.



**Easy to hold**

Head is mushroom-shaped and it is easy to carry.



**Easy to hang**

The holder designed with the best angle has achieved the ease of hooking pipette and secure holding.



**Smooth rotation**

Holder joint moves smoothly and revolves without frustration.



**Compact**

It is space-saving 6-holder carrousel stand. It is compact and stable.

\*The product package does not include pipettes.

Cat. No.	Item	Packaging
NT-STD	NEXTY Pipette Carrousel Stand	1 unit


# Repetitive Pipette

# Peggy

No electric power required at all!

Aspirate the sample all at once, dispense its necessary volume repetitively.


## Setting aspiration volume



Aspiration can be set in every 1 mL for setting range 1 mL ~13 mL.


Note: 13 mL can be dispensed in case that Watson 10 mL Disposable Pipette is used.

## Setting dispensing volume




Dispensing can be set in every 0.5 mL for setting range 0.5 mL ~5 mL.

## Flexible nozzle-angle



Adjustable nozzle angle in every 10 degrees up to 50 degrees.

## Simple maintenance




Easy to remove & set spare parts

- Nozzle cover: Autoclavable
- Joint packing: Autoclavable
- Water filter: Replaceable



Cat. No.	Item	Aspiration Volume max.	Aspiration Volume min.	Dispensing Volume max.	Dispensing Volume min.	Packaging
PG-01	PEGGY-01 Repetitive Pipette "Peggy"	13 mL	1 mL	5 mL	0.5 mL	1 unit/box



PG-01-JP1  
Joint packing

PG-01-JP2  
Joint packing  
(non filter type)

PG-01-SF1  
Replaceable  
water filter

## Consumables for Peggy


Joint packing is autoclavable.

Cat. No.	Item	Packaging
PG-01-JP1	Joint packing for Peggy	1 piece/bag
PG-01-JP2	Joint packing for Peggy (non filter type)	1 piece/bag
PG-01-SF1	Replaceable water filter for Peggy	5 pieces/bag





Repetitive Pipette "Peggy"



for  
**10 mL**  
5 mL  
pipette

**Precaution in case of using 5 mL disposable pipette:**  
Please DO NOT set aspiration volume exceeding 5 mL in case that 5 mL disposable pipette is used. It may cause an accidental liquid suction leading to malfunction on Peggy.

**Aspiration**

**1** Set the aspiration volume using the aspiration-volume-setting dial.  
*Note: DO NOT set volume exceeding the volume of disposable pipette.*

**2** Attach the Disposable Pipette to the Peggy nozzle.

**3** **Point** Hold the Peggy vertically. Push the aspiration button slowly to aspirate the sample. After aspiration, slowly lift the Peggy out of the test tube.

**Dispensing**

**1** Set the dispensing volume using the dispensing-volume-setting dial.

**2** **Point** Make sure that the dispensing button returns all the way to the top after dispensing. Push the dispensing button one by one. Hold down the dispensing button until the solution no longer comes out of the tip. Each time the button is pushed, the set volume is dispensed.

**Point** Hold the Peggy vertically.

Dispensing volume can be changed for each dispensing.

\*Specifications are subject to change for further improvement without prior notice.



Electric Multi Dispenser "Black Bird"

### Electric Multi Dispenser

# Black Bird

The dispensing stops exactly at the set volume!

Multi-functional Dispenser for:

- Aspiration/Dispensing,
- Repetitive dispensation,
- Continuous aspiration



**Aspiration/Dispensing speed is adjustable**  
**Gentle to tissues, Applicable to high viscosity solution**

Can set the aspiration/dispensing speed to prevent tissue damage or spattering.  
 Can aspirate a high-viscosity solution in low-speed mode.



**Can aspirate/dispense without visual check**

The aspiration and dispensing will automatically stop at the set volumes. (No need to check the scale of the disposable pipette.)  
 Once button is pressed, aspiration/dispensing starts automatically. This lessens your finger stress.



Aspiration/Dispensing Button

#### Black Bird tolerance range

0.5 mL	1 ~ 3mL	Over 3 mL
± 20%	± 10%	± 5%

\*You can set the volume from 0.2 to 0.5 mL, but this range is out of the tolerance guarantee range.  
 It is desirable that black bird is used as dispenser for volume of milliliter.  
 When dispensing of small volume, fine adjustment by using scale of disposable measuring pipette is needed.  
 When dispensing liquid quantitatively, it is recommended to aspirate more 10% extra volume than aspiration volume you want.  
 Example) If you wanted to dispense each 1 mL liquid at 10 times, you need aspiration of 11 mL liquids and dispensing each 1 mL liquids.

Cat. No.	Item	Aspiration Volume		Dispensing Volume		Packaging
		max.	min.	max.	min.	
BB-EMD1	Electric Multi Dispenser "Black Bird"	50 mL	0.2 mL	50 mL	0.2 mL	1 unit/box



\*Specifications are subject to change for further improvement without prior notice.



### Variable nozzle angle

The angle of the Black Bird nozzle is adjustable. Easy to operate on a clean bench with limited handling space, etc.



### Three-mode available

Black Bird has aspiration/dispensing mode, continuous dispensing mode, and continuous aspiration mode.

**Note:** Operator needs to change into the metal nozzle before changing mode from aspiration/dispensing to continuous dispensing or continuous aspiration.



### Easy to operate

The aspiration/dispensing volume of Black Bird can be easily set on the operation panel. Volumes can be set in every 0.2 mL. Frequently-used settings can be preset and reproduced.



### Double power supply drive

Either an AC adapter or a battery (AA cell) is available as power supply in Black Bird. (Battery visible alarm equipped.)



AC adapter [ Type A ]

An electrical current adapter

[input] AC 100V~240V/50Hz~60Hz 300mA  
[output] DC 5V 2A

B BF B3 C SE O

\*Please choose from these adapter at the time of an order.

### Parts for Black Bird

Cat. No.	Item	Packaging
BB-CB	Catch bottle	1 piece
BB-CAP	Catch bottle cap (with Silicon tube, outside diameter: 4 mm)	1 piece
BB-NA	Metallic Nozzle	1 piece
BB-JP	Joint packing for using measuring pipette	1 piece
BB-ST50	Silicon tube (outside diameter: 3 mm, length:50 cm)	1 piece
BB-PR	Tube pump (with PharMed®Tube, Tube adapter, mini connector)	1 piece