


what goes into pcr?

PCR has a lot of moving parts, each with a specific role and purpose. Understanding how these parts contribute to the PCR reaction and modulate each other is instrumental to successful PCR optimization.

Picking a Polymerase

- for routine PCR processes. It features relatively high thermostability and quick processing speed,¹ but possesses low replication fidelity
- Pfu: Originally isolated from Pyrococcus furiosus, Pfu is highly thermostable. *Pfu* also has proofreading capabilities, giving it extremely high replication fidelity (1-2 errors per million base pairs¹). However, *Pfu's* drawback lies in its slower processing speed, which is roughly
- **KOD:** A recombinant polymerase derived from *Thermococcus kodakarensis* KOD1 strain, thermostable KOD possesses the fidelity of *Pfu* with much higher processing speeds.⁴ However, KOD encounters difficulties with amplicons greater than 5 kb, resulting in lower yields.⁵



Identifying Contaminants

- **Extraction reagents:** Common reagents used during lysis and extraction such as phenol and SDS inhibit PCR activity by degrading
- Proteinases: Left over from tissue digestions, proteinases target
- Calcium: Calcium ions inhibit DNA polymerase activity by

DNA polymerases replicate DNA sequences during the elongation phase, so their function is pivotal to PCR success. Different polymerases have their own incorporation speeds, maximum amplicon lengths, and error rates.

DNA Polymerase

Buffers

А

В

С

D

E

F

Buffers ensure an environment that promotes DNA polymerase activity, with specific compositions varying depending on DNA polymerase selection. Tweaking buffer composition by adding chemicals or co-solvents (such as bovine serum albumin) may help when facing difficult samples such as GC-rich templates.



Template DNA

Template DNA quality and quantity will determine the overall amplification output. DNA fragmentation, excess starting material, and the presence of undesired DNA fractions can all cause non-specific amplification, while low template DNA concentrations result in limited output.



Primers

Primers guide DNA synthesis, ideally ensuring that only the sequence of interest is amplified. Primer binding sites should be unique; primer thermokinetic properties should roughly match that of the template DNA; and no complementary sequences within a primer or primer pair should exist.



dNTPs

Free deoxynucleoside triphosphate (dNTP) molecules are the building blocks for DNA amplification. However, excess dNTPs can reduce sequence replication fidelity, especially when using nonproofreading DNA polymerases.



Magnesium Ions


Magnesium ions are cofactors for DNA polymerase, enabling dNTP incorporation during elongation. Mg2+ ions also facilitate primer-DNA template complexing. Elevated Mg2+ concentrations increase DNA polymerase activity, but this can lead to non-specific PCR products.





(G

Template DNA quality and quantity will determine the overall amplification output. DNA fragmentation, excess starting material, and the presence of undesired DNA fractions such as genomic DNA can all cause nonspecific amplification, while low template DNA concentrations result in limited output.



References

- competing with magnesium ions for bindin g. Increasing magnesium this situation.
- EDTA: Since it inhibits DNAses, EDTA is often found in storage
- Chaotropic salts: Employed to facilitate DNA adsorption to

- F.C. Lawyer et al., "High-level expression, purification, and enzymatic characterization of full-length *Thermus aquaticus* DNA polymerase and a truncated form deficient in 5' to 3' exonuclease activity," *PCR Methods Appl*, 2(4):275-87, 1993.
 P. McInerney et al., Error rate comparison during polymerase chain reaction by DNA polymerase," *Mol Biol Int*, 2014:287430, 2014.
 S. Bustin, *The PCR revolution: Basic technologies and applications*. Cambridge: Cambridge University Press, 2010.
 M. Takagi et al., "Characterization of DNA polymerase from *Pyrococcus sp. strain* KOD1 and its application to PCR," *Appl Environ Microbiol*, 63(11):4504-10, 1997.
 M. Nichika et al., "Long and accurate PCR with a mixture of KOD DNA polymerase and its expoundese deficient mutant enzyme." J.
- 63(11):4504-10, 1997.
 M. Nishioka et al., "Long and accurate PCR with a mixture of KOD DNA polymerase and its exonuclease deficient mutant enzyme," J Biotechnol, 88(2):141-49, 2001.
 C. Schrader et al., "PCR inhibitors occurrence, properties and removal," J Appl Microbiol, 113(5):1014-26, 2012. 5.
- 6

© 2020 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved. Merck, the vibrant M, Sigma-Aldrich, GenElute and SmartLyse are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources Lit No. MK_AG6899EN

Sigma-Aldrich.